

REPORT DOCUMENTATION PAGE					Form Approved OMB No. 0704-0188	
<p>The public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing the burden, to Department of Defense, Washington Headquarters Services, Directorate for Information Operations and Reports (0704-0188), 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to any penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.</p> <p>PLEASE DO NOT RETURN YOUR FORM TO THE ABOVE ADDRESS.</p>						
1. REPORT DATE (DD-MM-YYYY) 04/11/2012		2. REPORT TYPE Quarterly Technical Report			3. DATES COVERED (From - To) 01/01/2012-03/31/2012	
4. TITLE AND SUBTITLE Human Neural Cell-Based Biosensor					5a. CONTRACT NUMBER N00014-11-C-0011, Mod. P00001	
					5b. GRANT NUMBER	
					5c. PROGRAM ELEMENT NUMBER	
					5d. PROJECT NUMBER	
6. AUTHOR(S) Stice, Steven L. Chilton, Jamie					5e. TASK NUMBER	
					5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) ArunA Biomedical, Inc. 425 River Road Athens, GA 30602					8. PERFORMING ORGANIZATION REPORT NUMBER QTR-11102010.4	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES) Director, Naval Research Lab Attn: Code 5596 4555 Overlook Avenue, SW Washington, D.C. 20375-5320					10. SPONSOR/MONITOR'S ACRONYM(S) DoD/ONR	
					11. SPONSOR/MONITOR'S REPORT NUMBER(S) ONR 5252.235-9714	
12. DISTRIBUTION/AVAILABILITY STATEMENT Approved for public release; distribution is unlimited.						
13. SUPPLEMENTARY NOTES						
14. ABSTRACT As a more biologically relevant model of human physiology, human neural progenitor cells have a strong potential as in vitro, cell-based biosensors for environmental toxins. The overall goal of this project is to develop a human neural cell based biosensor using Aruna's neural cell lines. Aruna's neural cell lines, hNP1™ Human Neural Progenitor Cells and hN2™ Differentiated Neuronal Cells, are derived from human embryonic stem cells (hESC), and more recently Aruna's neural cells lines are being derived from human induced pluripotent stem cells (hiPSC). This report describes progress in these major areas: (1) production of induced pluripotent stem cell derived neural progenitor cells, (2) directed differentiation of progenitors into dopaminergic neurons and astrocytes and (3) assay development of Aruna's neural cell lines as sensor elements for neurotoxicity.						
15. SUBJECT TERMS neurotoxicity, biosensor, neural cell						
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT UU	18. NUMBER OF PAGES 4	19a. NAME OF RESPONSIBLE PERSON Steven L. Stice	
a. REPORT Unclassified	b. ABSTRACT Unclassified	c. THIS PAGE Unclassified			19b. TELEPHONE NUMBER (Include area code) 706-583-0071	

Reset



425 River Road
Athens, GA 30605

**Quarterly Report
Human Neural Cell-Based Biosensor**

Date: April 11, 2012

Reporting Period: January 01, 2012-March 31, 2012

Prepared for:
Office of Naval Research (ONR)
Director, Naval Research Lab
Attn: Code 5596
4555 Overlook Avenue, SW
Washington, D.C. 20375-5320

Contract Number: N00014-11-C-0011, Amendment/Modification P00001

Submitted by:
Dr. Steven L. Stice, Principle Investigator
ArunA Biomedical, Inc.
425 River Road
Athens, GA 30602
Phone: 706-583-0071
Fax: 706-262-2821
Email: ssstice@arunabiomedical.com

Distribution Statement A:
Approved for public release; distribution is unlimited.

UNCLASSIFIED

Distribution of Quarterly Report

ADDRESSEE	DODAAC CODE	REPORT ENCLOSED	NUMBER OF COPIES	
			UNCLASSIFIED/ UNLIMITED	UNCLASSIFIED/ LIMIT ED AND CLASSIFIED
Program Officer: Dr. Laura Kienker ONR Code: 342 E-Mail: laura.kienker@navy.mil	N00014	Full technical report	1	1
Administrative Contracting Officer: Office of Naval Research ONR 0254: Russelle Dunson 875 North Randolph St. Arlington, VA 22203-1995 E-mail: russelle.dunson@navy.mil	S1103A	SF 298 only	1	1
Director, Naval Research Lab Attn: Code 5596 4555 Overlook Avenue, SW Washington, D.C. 20375-5320 E-mail: reports@library.nrl.navy.mil	N00173	Full technical report	1	1
Defense Technical Information Center 8725 John J. Kingman Road STE 0944 Ft. Belvoir, VA 22060-6218 E-mail: tr@dtic.mil	HJ4701	Full technical report	2	2

Summary

As a more biologically relevant model of human physiology, human neural progenitor cells have a strong potential as in vitro, cell-based biosensors for environmental toxins. The overall goal of this project is to develop a human neural cell based biosensor using ArunA's neural cell lines. ArunA's neural cell lines, hNP1™ Human Neural Progenitor Cells and hN2™ Differentiated Neuronal Cells, are derived from human embryonic stem cells (hESC), and more recently ArunA's neural cell lines are being derived from human induced pluripotent stem cells (hiPSC). This report describes progress in these major areas: (1) production of induced pluripotent stem cell derived neural progenitor cells, (2) directed differentiation of progenitors into dopaminergic neurons and astrocytes and (3) assay development of ArunA's neural cell lines as sensor elements for neurotoxicity.

(1) Production of induced pluripotent stem cells (iPSCs)

We have successfully been able to expand our new iPSC-derived neural progenitor cells and have completed our first production run. We are currently characterizing their karyotype and expression of neural progenitor cell phenotypic markers.

(2) Directed differentiation

Successful development continues on the directed differentiation of neural progenitor cells into dopaminergic neurons and astrocytes. Initial work to optimize differentiation protocols (media formulations, etc.) for hNP1™ derived dopaminergic neurons has been completed. Populations of hNP1™ cells were successfully differentiated into mature neuronal populations demonstrating positive protein expression of dopaminergic markers. Current work is now concentrating on reliably reproducing these results, as well as translating these culture methods to iPSC-derived neural progenitor cells. For astrocyte differentiation of hNP1™ cells, media conditions are being further refined and gene and protein expression is being more fully characterized at various timepoints. Preliminary work has also begun on translating astrocyte differentiation protocols to iPSC-derived neural progenitor cells. Work continues on the written manuscript of these findings on astrocyte differentiation for future publication.

(3) Assay development of hNP1™ and hN2™ cells as sensor elements for neurotoxicity

Progress this quarter on assay development has focused on 1) high content image analysis of neurogenesis and 2) high throughput automated patch clamp systems.

High content image analysis of neurogenesis: We have continued our development of high content image analysis of neurogenesis using an ImageXpress® high content imaging platform. Using this assay platform allows us to monitor neural progenitor cell proliferation, as well as differentiation in terms of neurite outgrowth in both our hNP1™ and hN2™ cells. This assay allowed for the detection of neurotoxicants as well as positive and negative modulators of neuronal differentiation, with output parameters including the number of neurites, length, branches, etc. per cell or per field. Results were recently presented in both oral presentation and poster format at the Society of Toxicology Annual Meeting, March 11-15, 2012, San Francisco, CA. In continued HCL studies, we are currently expanding our test set of known inducers of apoptosis and mediators of neuronal toxicity to market our hNP1™ and hN2™ cell lines with high content imaging platforms as a screening tool for drug discovery and toxicology studies.

High throughput automated patch clamp systems: Development has resumed on the large scale screening of compounds for effects on ion channel activity through direct electrophysical measurements using an IonWorks Barracuda™ system. Using this high throughput automated patch clamp system, the electrophysiological and pharmacological properties of endogenous ion channels in hN2™ cells can be characterized by measuring ligand and voltage-gated ion channels in 384 parallel recording sites. Such a platform technology can accelerate drug discovery, as well as evaluate neurotoxicity of compounds that are ion channel targets. Previous results were presented at the Society for Neuroscience Annual Meeting, November 12-16, 2011, Washington, D.C. Results from continued studies are expected to be presented in the upcoming annual meeting of the International Society for Stem Cell Research in June 2012.